

## Mass Spectrometry

# Immunoprecipitation (IP) protocol

### 1. Equilibrating the beads

*Except for section 3, all steps should be performed on ice or at 4°C.*

Protein A beads: Affi-Prep Protein A (Bio-Rad) beads have a binding capacity of 5-10 mg/ml rabbit IgG.

Protein G beads: GammaBind Plus Sepharose (GE Healthcare) beads have a binding capacity of 25 mg/ml mouse IgG.

In both cases, for initial binding tests we recommend the use of 1 mg of antibody per ml beads.

Wash beads 2× with 10 bed-volumes (bv) TBST (1×TBS + 0.04% Triton X-100) in 1.5 ml Eppendorf tubes: resuspend beads by flicking. Pellet beads by spinning for 10 sec in picofuge. Remove supernatant (SN) using suction line fitted with syringe and 27G needle.

### 2. Binding antibodies to beads

Add 10 bv TBST to the washed beads, then add purified antibody to give desired concentration. Mix beads with rotation for 1 h at RT or o/n at 4°C.

Wash beads 3× briefly with 10 bv TBST.

### 3. Cross-linking of antibodies to beads (*this section only: at room temperature!*)

Wash the beads 3× briefly with 10 bv 0.2 M Na-borate, pH 9.2.

(NB: the Na-borate solution forms crystals; microwave to re-dissolve, then cool to RT.)

Prepare a *fresh* solution of 20 mM DMP (i.e. 5.2 mg/ml) in 0.2 M Na-borate pH 9.2.

DMP = dimethyl pimelimidate dihydrochloride (Sigma D 8388. NB: DMP is highly moisture-sensitive; take DMP pot out of -20°C freezer 30 min before use, to equilibrate to RT).

Add 20 bv of this solution to the beads and mix at RT for **30 min: no longer !!!**

Stop cross-linking reaction by washing beads 2× 10 min with 10 bv 0.2 M Tris-HCl, pH 8.0.

Wash beads 2× briefly with 10 bv TBST to re-equilibrate.

**4. Pre-elution with glycine** (removes non-crosslinked antibodies)

Wash beads 2× briefly with 10 bv 0.1 M glycine, pH 2.0, then wash 3× briefly with 10 bv TBST to re-equilibrate. Ab-coupled beads can be stored at 4°C in TBST + 0.05% NaN<sub>3</sub>.

**5. Cell extract preparation**

Clarify whole-cell extracts by centrifugation in 1.5 ml epp(s) at 14,000 rpm, 4°C for 10 min (S14). Transfer the SN to a fresh tube and spin again. The SN from this spin is the clarified extract.

**6. Extract pre-clearing**

Transfer clarified extract to a 15 ml Falcon containing 1/10<sup>th</sup> volume protein-A/G beads, pre-loaded and cross-linked to non-specific IgG, and equilibrated with TBST. Incubate with rotation at 4°C for 30 min. Pellet beads and transfer SN (pre-cleared extract) to a fresh tube.

**7. Immunoprecipitation**

Transfer pre-cleared extract to a 15 ml Falcon tube containing specific antibody-coupled beads. Aim to saturate the beads with extract (e.g. 20 mg protein per 100 µl beads). Incubate with rotation at 4°C for 30 min. Pellet the beads, remove the SN (keep for analysis if desired).

**8. Washing the beads**

*(parameters in step 1 may vary according to Ab and complex, but washes indicated in step 2 and 3 have to be carried out!)*

***Step 1:***

**5 standard washes:** to pelleted beads, add 20-100 bv ice-cold IP-buffer plus 0.1% NP-40, 1×PIM, 0.1 mM PMSF. Resuspend beads by flicking, incubate with rotation at 4°C for 5 min. Pellet beads, aspirate as much SN as possible using needle.

**Optional: 3 higher-stringency washes:** as above, but with higher [salt] and [detergent], (e.g. 0.45 M NaCl), according to the desired stringency.

***Step 2:***

**10 no-detergent washes:** *(detergent must be removed: it interferes with MS/MS analysis!)*

Resuspend beads in 10 bv ice-cold wash-buffer. Spin, completely remove SN.

**For on bead digest, please stop here at this step. Both the elution and the tryptic digest will be performed by the protein/MS facility.**

**If you continue with a glycine elution, please continue with step 3.**

**Step 3:**

**2 washes to remove pH-buffering agents:** Wash beads twice briefly with 20 bv of 150 mM NaCl to remove pH-buffering agents. Resuspend beads in NaCl, transfer to 0.5 ml epp. Pellet beads and remove SN.

**9. Elution of immunoprecipitated protein from the beads ('glycine elution')**

*Under the hood:* add one or two bv of 0.1 M glycine (pH 2.0) to the beads, mix by flicking for 3 min. Pellet beads, transfer SN into a fresh Axygen low-retention 0.2 ml PCR tube, and repeat elution twice. Check for low pH of the combined eluates. Neutralize pH by adding 1 M Tris-HCl pH 8.0 (approx 15 ul for 100 ul eluate, check with pH paper).

Neutralise eluted beads by washing 3× briefly with 10 bv TBST. Transfer 1/10<sup>th</sup> vol eluate to fresh 0.5 ml Eppi, add 4× sample buffer for SDS-PAGE silver analysis. Use remaining eluate directly for in-solution digest, or freeze in liquid nitrogen and store at -80°C.

**10. IP bead recovery**

Add 10 bv 0.1 M glycine (pH 2.0) to the beads, mix by flicking for 1 min. Pellet beads and aspirate SN. Repeat once. Neutralise beads by washing 3× briefly with 10 bv TBST.

Store beads at 4°C in TBST + 0.05% NaN<sub>3</sub> for future use.

---

<b>IP-buffer:</b>	<b>stock conc</b>	<b>per 200ml</b>	<b>per 500ml</b>
20 mM Tris-HCl, pH 7.5	1 M	4 ml	10 ml
150 mM NaCl	5 M	6 ml	15 ml
10% glycerol	87%	23 ml	57.5 ml
2 mM EDTA	0.5 M	0.8 ml	2 ml

**PIM (Protease inhibitor mix):**

**Leupeptin, Pepstatin, Chymostatin**, each prepared at 10 mg/ml in DMSO. Dilute 1/1000 to generate 1x PIM in extract or wash buffer.

**Wash-buffer:**

20 mM Tris-HCl pH 7.5, 150 mM NaCl