

Mass Spectrometry

Gel Digestion Protocol

Gel digest - Trypsin

1. Cut the gel:

- Cut out bands you prefer and then cut each band in 2-3 mm pieces, put them in a 0.6 ml PCR tube from AXYGEM.

2. Wash the gel pieces:

- Use the shaker for each wash step at RT for 10 min. Remove the supernatant after each step.
- 200 μ l 100 mM ABC.
- 100 μ l ACN, 100 μ l 100 mM ABC.
- Repeat the two steps above.
- Remove the supernatant and add 100 μ l ACN to shrink the gel-pieces.

Note: if the bands are heavily coomassie-stained, you might wash them over night in 100 μ l ACN, 100 μ l ABC 100 mM on the shaker at RT.

3. Reduction:

- Cover the gel-pieces with 100 μ l of 1 mg/ml DTT in 100 mM ABC.
- Incubate at 57°C for 30 min.
- Remove supernatant.

4. Alkylation:

- Cover the gel-pieces with 100 μ l of 28 mM MMTS in 100 mM ABC.
- Incubate at RT in the dark for 30 min.
- Remove supernatant.

5. Wash the gel pieces:

- Perform the same wash steps which you did in step (2). Remove supernatant and shrink the gel-pieces as in step (2).
- Dry the gel pieces in the speed-vac for 5-7 min
- Close tubes and flick them carefully – dry gel pieces are loose and move around. If they don't, put tubes back in the speed vac for a few min

Note: It is now possible to store the gel-pieces at -80°C.

6. Trypsin digest:

Use only “trypsin gold” from “Promega” and our stock-solution ($c = 100 \text{ ng/}\mu\text{L}$). Dilute it with 100 mM ABC to $c = 12 \text{ ng/}\mu\text{L}$.

This is a 1:8 dilution step:

=> 160 μL aliquot: 20 μL stock + 140 μL ABC.

=> 80 μL aliquot: 10 μL stock + 70 μL ABC.

=> 40 μL aliquot: 5 μL stock + 35 μL ABC.

- Cover the gel pieces with 20 μL diluted trypsin/ABC solution and put it at 4-5°C for 5 min. If the gel-pieces suck in all of the solution, pipette another few μL of diluted trypsin/ABC solution on them and let them rest again for 5 min.
- Remove the supernatant and add the same volume of 100 mM ABC as you used before with the trypsin solution (20 μL or more).
- Incubate overnight at 37°C.

7. Extract gels:

- Pipette the supernatant (~20 μL) into a 0,2 ml PCR tube and store it at 5°C.
- Cover the gel pieces with 20 μL 5% formic-acid and sonicate them for 10 min in a *cooled* ultrasonic bath.
- transfer the 20 μL supernatant to the 0,2 ml PCR tube from the 1st step and put the tube in the fridge again.
- Once more, cover the gel-pieces with 20 μL 5% formic-acid and sonicate them in a *cooled* ultrasonic bath.
- Again, transfer the supernatant to the 0.2 ml PCR tube.
- Discard the gel pieces.

You should now have a 0.2 ml PCR tube with ~ 60 μL of your gel digest.

Gel digest - Chymotrypsin

The steps are the same as in the “trypsin digestion protocol”, there is only one important difference in step Nr. 6:

6. Chymotrypsin digest:

- Incubate at 25°C for 5h.

Gel digest - Subtilisin

1 – 5. Prepare the gel pieces:

- Prepare (cut, wash, reduce, alkylate, wash) the gel-pieces as described in the steps 1-5 of the “trypsin digestion protocol”.

6. Subtilisin digest:

Subtilisin has to be prepared freshly:

- Prepare the Urea/Tris dilution solvent by mixing 9 vol. parts of 6 M Urea with 1 vol. part of 1M Tris (eg.: 1800 μ L Urea + 200 μ L Tris).
- Prepare 5 mg subtilisin in 500 μ L 1mM HCl ($c = 10 \mu\text{g}/\mu\text{L}$).
- Dilute this 1:5 with 1mM HCl. ($c = 2 \mu\text{g}/\mu\text{L}$).
- Dilute this 1:160 with Urea/Tris dilution solvent ($c = 12.5 \text{ ng}/\mu\text{L}$).
- Dilute this 1:2 with Urea/Tris dilution solvent ($c = 6 \text{ ng}/\mu\text{L}$).
- Cover the gel-pieces with 30 μ L subtilisin solution and put it at 4-5°C for 5 min. If the gel pieces suck in all of the solution, pipette another few μ L of subtilisin solution on them and let them rest again for 5 min.
- Remove the supernatant and add the same volume of 100 mM ABC as you used before with the subtilisin solution (30 μ L or more).
- Incubate at 37°C for 1h on the shaker.

7. Extract gels:

- Pipette the supernatant (~30 μ L) into a 0,2 ml PCR tube and store it at 5°C.
- Cover the gel pieces with 20 μ L 5% formic acid and sonicate them for 10 min in a *cooled* ultrasonic bath.
- After sonication, transfer the 20 μ L supernatant to the 0,2 ml PCR tube from the 1st step and put the tube in the fridge again.
- Once more, cover the gel-pieces with 20 μ L 5% formic acid and sonicate them in a *cooled* ultrasonic bath.
- Again, transfer the supernatant to the 0.2 ml PCR tube.
- Discard the gel pieces.

You should now have a 0.2 ml PCR tube with ~ 70 μ L of your gel digest.

Gel digest - LysC

The steps and amounts are the same as in the “trypsin digestion protocol”